

Fecal sample into dry sample container

This instruction is followed for these tests:

- VirVi (2975)
- F-Elast1 (1718)
- F-VirNhO (14157)
- F-ParaNhO (13844)
- SienVi, feces (2631)

Each of these tests requires a separate sample container.

Materials

1. Plastic or cardboard container.
2. Container for fecal sample (Figur 1).
3. Plastic bag and sticker for name.

Procedure

1. Defecate into the plastic or cardboard container or on toilet paper.
2. Open the container for the fecal sample. Using the spatula attached to the top, take about one spoonful of feces into the container (do not fill the container).
3. If you have diarrhea, take one tablespoonful into the container.
4. Close the container tightly.
5. Please be careful so that the outside of the container does not get sample on it.
6. **Write your name, social security number, sampling date and sampling time on the sticker that came with the materials from the lab.**
7. **Attach the sticker to the sample container.**
8. Put the sample container into the plastic bag.
9. Keep the sample container cool, preferably at refrigerator temperature.



Figur 1

Take the sample container to the laboratory at the latest on the day after sampling.

Please contact the department or outpatient clinic if you have questions regarding the sampling procedure for this examination.

You may ask for the result of this examination at the hospital unit where you are being treated

Reviewed 17.7.2023 XX/PL